

Introduction to Single-Cell RNA-seq

The Data Lab

What can bulk RNA-seq vs single-cell RNA-seq help us determine?




Image from Tanaka *et al.* 2018 <https://doi.org/10.1038/s41551-017-0162-1>.

What can bulk RNA-seq vs single-cell RNA-seq help us determine?




Image from 10X Genomics blog: <https://www.10xgenomics.com/blog/single-cell-rna-seq-an-introductory-overview-and-tools-for-getting-started>

Single-cell RNA-seq quirks


Less starting material means:

- More PCR amplification (*and its associated biases*)
- More zero counts
 - Biology - Not every gene is expressed in every cell
 - Technical - Biased capture methods, Sequencing every RNA in every cell requires a lot more sequencing

Choi *et al.* (Genome Biology, 2020) <https://genomebiology.biomedcentral.com/articles/10.1186/s13059-020-02103-2>




Single-cell gene mean density graph




Single-cell gene mean density graph

A lot of genes that are low abundance



Single Cell Basic Set-ups

1. Tag-based scRNA-seq




Tag-based separation of cells' data *after* sequencing

Example: 10X Genomics Chromium

Zheng *et al.* 2017

<https://www.ncbi.nlm.nih.gov/pubmed/28091601>

2. Full-length scRNA-seq



Physical separation of cells *before* sequencing

Example: Smart-seq2

Picelli *et al.* 2014


<https://www.nature.com/articles/nprot.2014.006>

Cell Barcodes + Unique Molecular Identifiers (UMIs) are used to label individual transcripts

Each droplet contains 1 cell, all with the same cell barcode

Within each droplet, each bead contains millions of distinct UMIs


Each transcript within a cell is tagged with a cell barcode and unique molecular identifier (UMI)




Unique Molecular Identifiers (UMIs):

a 'snapshot' of the original molecules in the *pre-amplified* cell

1.




2.




Unique Molecular Identifiers (UMIs):

a 'snapshot' of the original molecules in the *pre-amplified* cell


1.



2.




3.




Unique Molecular Identifiers (UMIs): a 'snapshot' of the original molecules in the *pre-amplified* cell


1.




2.



3.




4.



Tag-Based scRNA-seq

Pros:

- Can profile up to millions of cells.
- Takes less computing power.
- File storage requirements are smaller.
- Much less expensive.




Tag-Based scRNA-seq

Pros:


- Can profile up to millions of cells.
- Takes less computing power.
- File storage requirements are smaller.
- Much less expensive.


Cons:

- More intense 3' bias because sequencing is not bidirectional.
- Coverage is generally not as deep as full-length scRNA-seq.





Single sample scRNA-seq overview





- Mapping & Quantification
 - Align reads and generate a gene by cell counts matrix
 - Cell Ranger, 10x supported tool, is popular mainly for being relatively user friendly, but is very slow (aligns to the entire genome)
 - Alevin is a faster salmon-based tool (aligns to the transcriptome)
- Import to R
 - Store gene by counts matrix as a SingleCellExperiment object
 - SingleCellExperiment objects are used to hold additional information about cells and genes in addition to the counts matrix



- Cell-level statistics
 - Calculate quality metrics for each cell
 - total UMI count
 - number of detected genes
 - mitochondrial content
- Filter cells
 - Remove low quality cells from counts matrix
- Normalize counts
 - Use size factors to minimize technical differences and maximize biological differences across cells
 - Log transform counts data

Some of the many resources for you in `00-scRNA-seq_introduction.Rmd`

- [Hemburg lab scRNA-seq training course](#)
- [ASAP: Automated Single-cell Analysis Pipeline is a web server that allows you to process scRNA-seq data.](#)
- [Smith. Unique Molecular Identifiers – the problem, the solution and the proof - article on background of UMIs](#)
- [Literature on technologies](#)